We bring quality to light.





CAS 140D IR: New High Resolution Infrared Array Spectrometer



Infrared Sensing Application - Overview

NIR SWIR

850nm

- Flood Illumination
- Eye-tracking

905nm

LiDAR

940nm

- Structured light
- Flood illuminator
- Proximity
- Eye-tracking
- LiDAR
- Heart rate
- SpO₂

1380nm

- Under display proximity sensors
- Fingerprint on display sensor
- World facing ToF

1550nm

LiDAR (fiber laser & EEL)

1465nm

Blood glucose

1050nm/1065nm

- Lactate
- Food analysis
- Body hydration
- Ethanol



F/W Facing



LiDAR



Gesture-tracking



DMS



Eye-tracking




Under-display



Bio-sensing



Calibration Chain





Metrological Excellence of Instrument Systems


Calibration competence

- Traceable calibration & verification of instruments and reference sources
- Device optimized audit concepts
- ISO 17025 accredited device testing
- Metrology lab in Shanghai (mirror site)
- Service network in Asia (China, Korea, Taiwan, Vietnam) – part of IS quality management and traceability network

Audit

- ACS sources in factory or lab (identical to IS in-house OQC sources)
- Self-absorption correction for different DUTs
- Customer support for implementation of audit processes

Ensured traceability and error budget






ACS IR Calibration and Audit Standards



ACS 570: Advanced calibration standards

- Principle: Highly stable infrared source based on LED technology
- Wavelength range: λ_{peak} 860nm, 950nm, 1200nm, 1300nm, 1380nm, 1450nm
- Compact design for inline applications
- Recalibration interval: 100 h
- Scope: Reference value for radiant flux in infrared
- Highlight: Automatic detection of driving conditions with ACU-100 control unit
- Highlight: Specific audit sources for spectrometer, photodiode and 2D cameras





Technical Key-Features CAS140D IR

- ✓ High-End InGaAS Sensor 512 Pixel, ultra high sensitivity, low noise and dark current
- ✓ Thermal Stabilization with TEC @ -10°C Low dark current, perfect long term stability
- ✓ Faster electronics, faster processing of measuring data Reduced measuring times, increased productivity and throughput
- ✓ Improved Spectrograph Design Minimum straylight and enhanced throughput reduction of optical artifacts like reflections, increased repeatability







Comparison	CAS 140CT	CAS 140D	
Min. integration time	10ms	1ms	
Min. scan time	16ms	9ms	1



Technical Key-Features CAS140D IR

Automatic accessory recognition
Failsafe probe recognition by optical PLG connector (plug & play)

Error and fail safety

Onboard memory. Calibration files (*.isc / *.ini) are stored in the CAS 140D with automated software check

Enhanced dynamic range
8-position filter wheel with option for up to 8 OD filters

Mechanical robustness

- Magnetic lockable dust cap for PLG connector
- Robust sheet metal housing
- Exchangeable dust filters
- Vibrational isolation / damping of the internal optics

Three color status bar








CAS 140D IR: Higher Sensitivity








CAS 140D IR: Lower Stray Light







Summary

- ✓ Instrument Systems provides high precision spectrometers, systems and software solutions for spectral light measurement, with their latest addition being the CAS 140D IR.
- ▲ All spectrometers and other products are calibrated with traceability back to national standards. Instrument Systems' lab is ISO 17025 accredited.
- Customers get close support with their audit concept in order to monitor repeatability and device drift. IS sells audit sources which closely match the customer DuTs.



THANK YOU for your attention!