

“Fiber optic sensing concepts for Minimally Invasive Surgery (MIS)”

June 12, 2024

Dr. Christian Voigtländer

R&D Manager
FBGS



Outline

- **About FBGS**
- **Minimally Invasive Surgery (MIS)**
- **Catheter navigation**
 - 3D Shape Sensing of catheters
 - Navigation using pre-operative CT-imaging data
 - Absolute shape tracking using augmented reality

Outline

- **About FBGS**
- Minimally Invasive Surgery (MIS)
- Catheter navigation
 - 3D Shape Sensing of catheters
 - Navigation using pre-operative CT-imaging data
 - Absolute shape tracking using augmented reality

we are

a developer and manufacturer of fiber optic sensor components and fiber optic sensor solutions



We are vertical from raw materials to software



we do

Temperature sensing



Shape sensing



Strain sensing



Force sensing



Pressure sensing



for

Medical



Process Industry



Civil engineering



Transportation



Energy



Expert in tailored fiber optic sensors and interrogation devices


High level of applications know-how
Customer oriented

we speak

18 languages and our team originates from 20 countries



Locations



FBGS North America (Montreal)
Sales (North American market)

FBGS Belgium (Geel)
Sales and applications
Development and assembly of measurement systems (interrogators)
Sensing Solution Engineering and R&D Work



FBGS Germany (Jena)
Company Headquarter
DTG® and FSG® production
R&D location for specialty fibers and FBGs




FBGS China (Suzhou)
Sales and applications (Chinese customers)




Fiber Bragg Gratings

Permanent periodic modulation of the refractive index in the fiber core:



BRAGG CONDITION: $\lambda_B = 2 n \Lambda$



$$\lambda_B (T, \varepsilon) = 2 n(T, \varepsilon) \Lambda(T, \varepsilon)$$

FBGS base technologies



Draw Tower Grating® technology

FemtoSecond Grating® technology

FBGscan interrogator platform



preform diameter
fibre diameter
(typically 125 µm)
drawing velocity



Features


- Stand-alone FBG measurement system
- Up to 40 sensors per optical channel
- No need for time of flight corrections to maintain wavelength accuracy
- Depolarized light source to reduce birefringence induced noise effects
- High dynamic range $\geq 30\text{dB}$
- Simultaneously detection of sensors belonging to the same optical channel
- Built-in processing board for calculating engineered data

FBGS product line and technologies


FBG sensor components

Interrogators

Software



Draw Tower Gratings (DTG®)



FemtoSecond Gratings (FSG®)



All Grating Fibers (AGF®)



Sensors



Measurement devices




Software


- High-strength FBG sensors produced during the drawing of an optical fiber



- High-strength FBG sensors produced by writing through the fiber coating



- Optical fiber containing densely spaced low reflective DTG®s ($R < 0.1\%$)
- More than 95% of the fiber contains gratings



- Variety of high-quality FBG sensors available for strain, force, temperature, pressure and shape measurements



- High-precision measurement devices for all FBG sensors
- Available up to 15 optical channels
- For strain, force, temperature, pressure and shape



- ILLumiSense software (embedded SW)
- Software for, data acquisition, processing and visualization
- Rest API
- Linux support



Key advantages fiber optic sensors



RELIABILITY

- Passive component
- Long lifetime (>20 years)
- No corrosion
- Stable over time (No calibration required)
- Cables and connectors are telecom grade

IMMUNITY

- **Immune to electro-magnetic radiation & radio frequency interference**
- Immune to high voltage discharge
- Explosion safe

SIZE

- Fiber is also the sensor
- **Lightweight & small diameter (< 1/4 mm)**
- **High integration and embedding capabilities**

PERFORMANCE

- **10's, 100's or 1000's of sensors in 1 fiber**
- **Less cables**
- Easy installation
- High fatigue resistance
- Long distance measurements (20+ km)

OMNI-FUNCTIONALITY

- Measurement of
 - Strain
 - **Temperature**
 - Vibration
 - Acoustics
 - **Pressure**
 - **Force**
 - **Shape**
 - Acceleration
 - Displacement
 - **Flow**
 - Liquid level

Outline


- About FBGS
- **Minimally Invasive Surgery (MIS)**
- Catheter navigation
 - 3D Shape Sensing of catheters
 - Navigation using pre-operative CT-imaging data
 - Absolute shape tracking using augmented reality

Minimally Invasive Surgery (MIS)


Minimally invasive surgery refers to **surgical procedures performed with minimal damage to the body's tissues**. These procedures are designed to minimize the size and number of incisions, reduce trauma to surrounding tissues, and promote faster recovery and reduced postoperative pain for patients

Example: Atrial Fibrillation (AFib)

- Replacement of **regular** heart rhythm with an **irregular fibrillatory** heart rhythm
- Possible treatment:
 - **Catheter ablation** (interventional procedure)
- **Regions responsible for the irregularly generated electrical impulses are ablated to electrically insulate them**



Courtesy of the MAYO foundation for medical research and the U.S. National Library of Medicine



Current challenges, limitations and opportunities

① Patients

- ❖ Overall instrument interaction
 - Vessel walls
 - Soft tissue
- ❖ Ablation tip contact forces
 - Perforation
 - Esophageal injury

② Interventionalists

- ❖ Ionizing radiation (fluoroscopy)
 - Cancer
 - Orthopedic injuries
- ❖ Limited visualization
 - 2D grey-scaled images

③ Overall procedure

- ❖ Safety and precision
- ❖ Sensing & visualization

Opportunities of FO sensing in MIS



- Provide 3D shape information of catheter(s)
 - Assist with localisation/navigation of catheter inside body
 - Measure intended or non-intended interaction forces between catheter and body
- } Reduce exposure to ionizing radiation

Outline

- About FBGS
- Minimally Invasive Surgery (MIS)
- **Catheter navigation**
 - 3D Shape Sensing of catheters
 - Navigation using pre-operative CT-imaging data
 - Absolute shape tracking using augmented reality



Shape sensing using MCF-DTG®



Differential Frenet-Serret formulae


$$\frac{dT}{ds} = \kappa N,$$

$$\frac{dN}{ds} = -\kappa T + \tau B,$$


$$\frac{dB}{ds} = -\tau N.$$

With $\tau = \frac{d\theta_b}{l_z}$,


Measure wavelengths




Compute strain



Interpolate strain



Compute \kappa and \theta_b



Compute 3D shape



Reference: Al-Ahmad, O., Ourak, M., Van Roosbroeck, J., Vlekken, J., Vander Poorten, E. (2020). Improved FBG-Based Shape Sensing Methods for Vascular Catheterization Treatment. IEEE Robotics and Automation Letters, 5 (3), 4687-4694. doi: 10.1109/LRA.2020.3003291



Shape sensing for endoscope navigation

Endoscope navigation

- Navigation of endoscope in the body
- Avoidance of endoscope looping
- Insertion length measurement
- Tip torture detection




Outline

- **About FBGS**
- **Minimally Invasive Surgery (MIS)**
- **Catheter navigation**
 - 3D Shape Sensing of catheters
 - Navigation using pre-operative CT-imaging data
 - Absolute shape tracking using augmented reality

Navigation using pre-operative CT-imaging data




- **Computed Tomography scan (CT-scan)**, is a medical imaging technique that uses X-ray technology and advanced computer processing to create detailed 3D images of the body, including the bones, muscles, organs and blood vessels
- Segmentation technologies can be used to extract 3D shape information of e.g. blood vessel



Example of Coronary segmentation from CT-scan image
(Courtesy: Adriaan Coenen, Medical Center, Rotterdam)

- CT scans are frequently used as part of preparation for surgery or treatment.
- **3D shape information from CT-scan can be used as a reference shape to localize the catheter**

Navigation using pre-operative CT-imaging data





Patent application: Determining position or insertion length of an elongated device (US 18/164,817 and EP20707122.6).

Outline


- **About FBGS**
- **Minimally Invasive Surgery (MIS)**
- **Catheter navigation**
 - 3D Shape Sensing of catheters
 - Navigation using pre-operative CT-imaging data
 - Absolute shape tracking using augmented reality

Absolute shape tracking using augmented reality



- **The depth camera allows to perform spatial mapping**
= the ability to create a **3D map of the environment** and is an inherent part of augmented reality technology.
- Spatial mapping makes it possible to place objects on real surfaces.
- Relative position in 3D space calculated from position markers robe
- **Hologram of the shape probe is projected in space matching the physical presence of the shape probe.**

Demonstration absolute shape tracking using augmented reality



It's a very promising technology because on the one hand, it allows you to visualize the catheter in 3D in relation to a patient.

Questions?

Dr. Christian Voigtländer (R&D Manager)

- cvoigtlaender@fbgs.com
- www.fbgs.com

- **FBGS International NV**
Bell-Telephonaan 2H
2440 Geel (Belgium)
- **FBGS Technologies GmbH**
Franz Loewen Str. 3
07745 Jena (Germany)

