

Hyperspectral Imaging for Pharma Applications

Julio Hernández

Hyperspectral Applications Manager julio@neo.no hyspex.com

Norsk Elektro Optikk AS

About Norsk Elektro Optikk (NEO): with over 35 years since its establishment, NEO is the largest private R&D company in the field of electro-optics in Norway. Together with its subsidiaries in Norway and Sweden, NEOs activities focus on the development of state-of-the-art products for the international market. HySpex is the Hyperspectral Imaging Division of NEO.



HySpex is a research-based manufacturer of hyperspectral cameras with over 15 years of experience providing solutions for scientific and industrial applications. HySpex sensors are recognized for their scientific-grade accuracy and are used worldwide in the most demanding applications in the VNIR-SWIR spectral range.

HySpex *Classic* - Scientific



HySpex Mjolnir - UAV



HySpex Baldur - Industrial






How can we contribute to improve pharma processes?


Hyperspectral imaging

- Main advantage: identifies spectral properties distributed on an area with high spatial resolution in real time.
- Classification: pharmaceutical components vs foreign objects and/or contaminants.
- Identification: specific API or mixture components; distribution on samples.
- Quantification: actual API content and concentration on samples.

HySpex

- Flexible instrumentation: multiplatform deployment.
- Feasibility studies: is HSI suitable for your process?
- Turn-key solutions: custom development.








Example: API identification and distribution in mixes




Sample preparation



Imaging



Spectral library or sample characterization



PLS modeling



Our technology

- The HySpex hyperspectral cameras are characterized by high spatial resolution, high spectral fidelity and low optical aberrations.
- High spectral resolution cameras in the VNIR-SWIR spectral regions (400-2500 nn) with down to micrometer spatial resolution.
- Our developments for high-end, scientific applications has allowed an efficient technology transfer to industrial systems with high repeatability and accuracy.
- High speed sensors and SDK allow for a smooth real-time implementation on a variety of research and industrial environments.
- Robust optomechanical design allows for repeatability and scalability on a wide range of applications





Thank you!



Questions and Contact hyspex@neo.no www.hyspex.com

