


EPIC MEETING PRESENTATION

Luca Cerri

V.P. Marketing and Scientific Director Dermatology and Aesthetic Medicine



Quanta System

Founded in 1985

Active in Dermatological, Surgical and Art Restoration




New, state-of-the-art 13,000m² facility

ISO 13485:2003; 13485:2012; 13485:2016; AC:2012; 9001:2008

EC Certificate: Full Quality Assurance System, Directive 93/42/EEC

EC Certificate: Production Quality Assurance System, Directive 93/42/EEC


Over 230 employees supporting 110 distributor partners
globally committed to world-class support








Medical Portfolio

Market leader in medical field thanks to outstanding clinical performances based on exclusive solution





Mixed technology

Mixed Technology working scheme




Hair Removal






Mixed technology

Ablative and non ablative lasers




Emission of 10600 nm only

*Standard
emission*




Emission of 1540 nm only





SEQUENTIAL emission of 10600 nm and 1540 nm with independent pulse width, power and variable delay between the two pulses



*MT
modality*



SIMULTANEOUS emission of 10600 nm and 1540 nm with independent pulse width and power



Ablative and non ablative lasers




**Fractional Resurfacing
Mixed 1540+10600 nm**



Before




After 60 days



Discovery Pico – Fractional Effect




Quanta System




Acne Scars



Fractional picosecond laser treatment, 1 session





THANK YOU

